





BumbleBee™ SAMPLE HANDLER

SPEED

Fast transfer rates as low as 3 seconds per sample, ideal for plate reformatting and hit picking

VERSATILE

Systems are configurable from 1 to 16 single-channel pipettors and multi-channel pipetting heads. An example configuration includes 2 tube pickers, 6 single-channel pipettors, 1 96-well multi-channel pipetting head, and 1 384-well multi-channel pipetting head.





General Specifications

Single-Channel 70uL, 250uL, and tube picking options

Multi-Channel Head Up to 2, compatible with Agilent heads

Software VWorks Automation Control Software

Serviceability <60 second pipette/syringe module swaps