



HiG™ AUTOMATED CENTRIFUGE

CAPABILITY

Up to 100 grams of imbalance tolerance, 350 gram payload, 5000 RCF, and active cooling to maintain ambient conditions.

FLEXIBILITY

Low vibration, easy to use standalone interface, and a .NET API for integrations

RELIABILITY


Proven reliability with minimal downtime during routine PM's and service



www.bionexsolutions.com

©2024 BioNex Solutions Inc.

LIT-50005 Rev G Effective 2024-05



General Specifications

Imbalance	100g for HiG3 and 50g for HiG4
Active cooling	Maintains ambient conditions (HiG4)
Labware	SBS up to 61mm tall
Electrical	100-240 VAC, 50/60 Hz, 15 A