



# VWORKS™ AUTOMATION CONTROL SOFTWARE

## EFFICIENCY


Native and direct controls for BioNex and Agilent devices and systems allows for fast and efficient workflow scheduling


## CAPABILITY

Custom driver development for new automation equipment is easily achieved with our rapid on-site development team or included SDK package

## FLEXIBILITY

True device pooling and a robust automated error recovery increases your walkaway time





### General Specifications and Features

<b>Computer Req</b>	Windows 10, 8 GB memory, 250 GB HD
<b>Scheduler</b>	Drag and Drop GUI for workflow design
<b>Packages</b>	SDK, Hooks, and Watcher
<b>Controls</b>	User management and robust logging